

Service Information

Singulex offers antibody derivatization services for you to conduct your own assay development projects. We adhere to the utmost quality standards in manufacturing, providing solutions designed specifically to answer your research needs.

Our Process

Send us your antibodies of choice and we will send you back capture antibody-coated microparticles and AlexaFluor®-labeled detection antibodies. Materials to be labeled can either be supplied by you or purchased by Singulex from an external vendor. Our derivatized reagents are optimized for use on the Erenna® Immunoassay System.


For best results, materials supplied are to be provided in carrier- and label-free buffer.

Minimum Capture Antibodies Required	Minimum Detection Antibodies Required	Derivatized Kit Equivalents
1 mg	300 µg	30-35 kits

Our Advantage

Singulex derivatization service is ideal for proof-of-concept assay development projects. Not only will your proprietary information be kept confidential, you will have access to quality reagents to conduct feasibility studies on your Erenna® Immunoassay System.


For additional assay development support beyond basic antibody derivatization, Singulex provides custom contract assay development and verification services to help you fast-track your biomarker research.



Our Promise

- quality labeling
- full confidentiality
- fast turn-around time
- expert phone support

We ship within 2-4 weeks after receiving your antibody and order in most cases.



Related Products and Services

Contact us today for related reagents and services to help you jump-start your assay development project!

Part Number	Product Name (Unit Amount)
CUSTOM0003	Derivatization for Coated Beads (up to 1 mg Ab)
CUSTOM0007	Derivatization for Detection Antibody (up to 1 mg Ab)
96-0005-00	Sample Testing Service
96-0003-00	Custom Assay Development Service
02-0001-05	10X Wash Buffer (1 L)
02-0046-02	Erenna® 20 Set Calibrators (1 box)
02-0111-00	Erenna® 10X System Buffer (1 L)
02-0063-02	Assay Buffer (1 L)
02-0297-00	Elution Buffer B (100 mL)
02-0368-00	Buffer D (100 mL)

Note: Singulex provides generic assay buffer and elution buffer for assay development. Purchaser must determine the suitability of the product(s) for their particular use.